

# Today's Agenda

- ✓ Look at two *M. tuberculosis* datasets
  - ✓ Run fastqc to look at fastq data quality
  - ✓ Trim poor quality reads with Trim Galore
  - ✓ Mapping module
-

# ILLUMINA DATA QC




# Illumina




fastq

# Nanopore



fast5

# Ion Torrent



bam



# Fastq format

```
1  @SEQ_ID
2  GATTGGGGTTCAAAGCAGTATCGATCAAATAGTAAATCCATTTGTTCAACTCACAGT
3  TT +
4  !"*((((***+))%%%++)(%%%).1***-+*"))**55CCF>>>>>CCCCCCC65
```

**Line 1** begins with a '@' character and is followed by a sequence identifier and an optional description (like a FASTA title line).

**Line 2** is the raw sequence letters.

**Line 3** begins with a '+' character and is optionally followed by the same sequence identifier (and any description) again.

**Line 4** encodes the quality values for the sequence in Line 2, and must contain the same number of symbols as letters in the sequence.

# Fastq format

## fastq header format (version > 1.8)

Sequence Header							+Sequence ID				
a	b	c	d	e	f	g	h	i	j	k	
@	HWI-ST486	:166:	C06K9ACXX	:7:	1101:	1443:	1995	1:	N:	0:	ACAGTG

**a. unique instrument name**

b. run id

c. flowcell id

d. flowcell lane

e. tile number within the flowcell lane

f. x-coordinate of the cluster within the tile

g. y-coordinate of the cluster within the tile

**h. the member of a pair, 1 or 2 (paired-end or mate-pair reads only)**

i. Y if the read fails filter (read is bad), N otherwise

j. 0 when no control bits are on

k. index sequence



[illegible]

Each sequence read is represented by 4 lines

ader

or

[illegible]

# Quality score interpretation

$$Q = -10 \log_{10} P \quad \longrightarrow \quad P = 10^{\frac{-Q}{10}}$$

Phred Quality Score	Probability of incorrect base call	Base call accuracy
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1000	99.9%
40	1 in 10000	99.99%
50	1 in 100000	99.999%

The quality (Q), also called Phred score, is the probability (P) that the corresponding basecall is incorrect.



# fast5 format

Binary file (not human readable)

Contains:

- Sequence of a read
- Raw signal data from pore
- Additional log files

Typically convert fast5 to fastq for downstream analyses





# BAM format for read data


**File format: SAM / BAM** (each line: one aligned sequence read)

The SAM/BAM file format is very powerful. It is unlikely that you will need to work with the contents of a SAM/BAM file directly, but it is very informative to visualize it in a viewer and it is a great format to do further analysis with. The format specifications are at <http://samtools.sourceforge.net/SAM1.pdf>. Below is a brief overview of the information contained in such files.

Binary Alignment Map format


Binary conversion of the Sequence Alignment Map (SAM) file

Typically convert bam to fastq for downstream analyses





# The “Golden” Rule



Having good quality fastq data is important!!










# FastQC: Per base sequence quality

## Good data

## Bad data



# MultiQC





# Quality score interpretation

$$Q = -10 \log_{10} P \quad \longrightarrow \quad P = 10^{\frac{-Q}{10}}$$

Phred Quality Score	Probability of incorrect base call	Base call accuracy
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1000	99.9%
40	1 in 10000	99.99%
50	1 in 100000	99.999%

The quality (Q), also called Phred score, is the probability (P) that the corresponding basecall is incorrect.

# Many tools for trimming

- Trimmomatic
- sickle
- fastP
- bbduk
- cutadapt
- Trim Galore

**Questions?**